



Alaris® Enteral syringe pump

Designed to deliver nutrition efficiently, the Alaris® Enteral syringe pump includes a large quick-to-read screen, intuitive keyboard, alarm beacon and familiar user interface. The Alaris® Enteral syringe pump recognises a wide range of enteral syringes and is clearly identified with its purple labelling for fast identification of enteral use.

Alaris® Enteral syringe pump



Designed to deliver nutrition efficiently, the Alaris® Enteral syringe pump includes a large quick-to-read screen, intuitive keyboard, alarm beacon and familiar user interface. The Alaris® Enteral syringe pump recognises a wide range of enteral syringes and is clearly identified with its purple labelling for fast identification of enteral use.

Alaris® Enteral syringe pump Specifications

Rate Range

5ml syringe	–	0.1 to 150ml/h
10ml syringe	–	0.1 to 200ml/h
20ml syringe	–	0.1 to 200ml/h
30ml syringe	–	0.1 to 200ml/h
50ml syringe	–	0.1 to 200ml/h

Volume Delivered

0.0ml - 9990ml

Volume To Be Delivered range

0.1ml - 1000ml, 1 min - 24hrs

Bolus Rate

5ml syringe	–	10 to 150ml/h
10ml syringe	–	10 to 200ml/h
20ml syringe	–	10 to 200ml/h
30ml syringe	–	10 to 200ml/h
50ml syringe	–	10 to 200ml/h

Max Purge Rate

500ml/h (limited to size of syringe)

Max Purge Volume Range

Purge volume range is 2ml

Volume to be delivered over time

Pre-defined volume can be set to deliver over pre-defined time.

Battery Type

Rechargeable sealed NiMH. Automatically charges when the pump is connected to AC power.

Battery Life

Mean Time To Battery Empty from fully charged at 5ml/h & 23°C under normal conditions is 6 hours.
(95% lower confidence interval of 5 hours 50 minutes)

Recharge Time

Charging takes 2½ hours from discharge to 90% charge.

Pumping Pressure Limits

11 Levels L0-L10 (50mmHg - 1000mmHg)

Alarms, Warnings and Prompts

Audible alarm and display prompts with comprehensive protection such as:

- Drive disengaged
- Check syringe
- Near end of delivery (NEOD)
- VTBD Done
- Internal malfunction
- Occlusion
- Battery low
- End of Delivery (EOD)
- Attention (nurse callback)

AC Power Requirements

115-230 VAC, 50/60Hz. 20VA (nominal)

Dimensions

310mm(w) x 121mm(h) x 200mm(d)

Weight

2.7 Kg

	5ml	10ml	20ml	30ml	50ml	60ml
CareFusion Enteral	0000ME00870	0000ME00871	0000ME00872	0000ME00873	✘	0000ME00874
Pentaferte Enteral	2022590	2022690	2022790	✘	✘	2022990
Vygon Enteral	✘	1015.102	1015.212	✘	✘	1015.602
Vygon A-VY Enteral	✘	✘	1015.213	✘	✘	1015.603
Medicina B.Tip	✘	✘	✘	✘	✘	Bladder Tip PE60B
Medicina Enteral	PE05	PE10	PE20	✘	✘	PE60
Nutricair	✘	SE10	SE20	✘	SE50	✘
Terumo Enteral	✘	✘	✘	SS+30EO	SS+50EO/C	✘

CareFusion International Headquarters
A-One Business Centre
Zone d'activités Vers-la-Pièce no 10
1180 Rolle
Switzerland
+41 21 556 30 00 tel
+41 21 556 30 99 fax

CareFusion South Africa
Unit 2 Oude Molen Business Park
Oude Molen Road
Ndabeni
Cape Town 7405
South Africa
+27 21 510 7562 tel
+27 21 510 7567 fax

CareFusion UK Ltd
The Crescent, Jays Close, Basingstoke
Hampshire RG22 4BS
United Kingdom
0800 917 8776 freephone
+44 (0)1256 330860 fax

CareFusion
PO Box 5527
Dubai
United Arab Emirates
Tel: +9714 2822 842
Fax: +9714 2822 914

© 2010 CareFusion Corporation or one of its subsidiaries. All rights reserved.
Alaris is a registered trademark of CareFusion Corporation or one of its subsidiaries. 0000CF00567 Issue 2